
HMF-IVDD-430

Manual de usuario de OncoAct

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2 Etiqueta

	V5.33-1.0		
	In Vitro Diagnostic medical device		Hartwig Medical Foundation Science Park 408 1098 XH Amsterdam www.Hartwigmedicalfoundation.nl
	(01)08720299486058(8012)v5.33-1.0		0344
	Instructions for use are supplied in electronic form instead of paper form. URL: https://oncoact.nl/manual/ Email: diagnosticssupport@hartwigmedicalfoundation.nl Device with internet access, web browser and PDF reader required for reading the manual. Paper instructions for use can be requested at no additional cost by contacting Hartwig using the indicated e-mail address and will be delivered within 7 days		

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An OncoAct report can be identified by the:

- Hartwig Medical Foundation logo at the top left corner on all pages of the report
- Title "Hartwig Medical OncoAct" at the top center of all pages of the report
- Signature of the Director Hartwig Medical Foundation on the last page of the report

One page summary with the most important results of the whole genome sequencing (WGS) analysis.

Primary tumor location and type as provided by the requesting medical expert.

Concise textual summary of the most relevant findings and their potential treatment options.

Overview of the main genomic tumor characteristics:

- Molecular tumor cell purity as measured using the sequencing data
- Molecular tissue of origin prediction
- Mutational burden status (low or high)
- Microsatellite status (stable - MSS, or unstable - MSI)
- Homologous recombination (HR) status (proficient or deficient)
- Tumor-associated viruses

More details are provided on page 5 and 6.

Overview of the main genomic tumor alterations:

- Genes with driver mutation(s)
- Genes with substantial copy gain (amplification)
- Genes that are completely lost in the tumor
- Genes that are completely disrupted in the tumor
- Gene fusions (in-frame and potential activating)

More details are provided on page 3.

Hartwig Medical OncoAct

OncoAct tumor WGS report

Summary

PRIMARY TUMOR LOCATION: Skin

PRIMARY TUMOR TYPE: Melanoma

The information regarding the primary tumor location and type, and the information related to the biopsy, is based on information received from the originating hospital.

Summary of most relevant findings

- Molecular tissue of origin prediction: Melanoma (likelihood: 99.6%)
- BRAF p.V600E activating mutation, possible indication for BRAF and MEK inhibitors
- CDKN2A p.A486T, p.G486S inactivation
- PTEN (negative) (0 loss)
- TERT (G-125, -124delC) promoter mutation
- TMB (13.7) positive, possible indication for checkpoint inhibitors

An overview of all detected cancer associated DNA alterations can be found in the report

Further interpretation of these results within the patient's clinical context is required by a clinician with support of a molecular tumor board

Tumor characteristics

Tumor purity: 99%

Molecular tissue of origin prediction: Melanoma (99.6%)

Tumor mutational burden status: High (13.7)

Microsatellite status: MSS (0.1)

HR Status: Proficient (0.00)

Virus: None

Genomic alterations in cancer genes

Genes with driver mutation: BRAF, CDKN2A, TERT

Amplified gene(s): None

Deleted gene(s): PTEN

Homozygously disrupted genes: None

Gene fusions: None

Pharmacogenetics

GENE FUNCTION

DPYD Normal Function

UGT1A1 NA

HLA Alleles

GENE GERMLINE ALLELE

HLA-A A*01:01

HLA-B B*40:02 [B*08:01

HLA-C C*07:01 [C*03:04

Germline results

Underlying data concerning cancer predisposition genes may be requested by a licensed clinical genetics laboratory after the patient has given informed consent.

Patient and sample details as provided by the requesting medical expert.

The status of the patient's genes involved in drug metabolism (pharmacogenetics) and status of the HLA alleles are summarized here. More details are provided on page 4.

This report is focused on the identification of all oncogenic driver alterations and the potential targets for therapy. Only tumor-associated genomic alterations in cancer associated genes are reported. The complete list of genes analyzed in this report can be found in <https://www.oncoact.nl/specsheetOncoActWGS>.

Germline results may be requested by a clinical geneticist, when informed consent was given by the patient.

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Genomic based therapy approaches

High level evidence

DRUG TYPE	TUMOR TYPE SPECIFIC	MATCH	LEVEL	RESPONSE	GENOMIC EVENT
BRAP inhibitor	Yes	Codon -600, Hotspot	A	▲	BRAP p.V600E
BRAP inhibitor, MEK inhibitor (Pani, MEK1 inhibitor, MEK2 inhibitor)	Yes	Codon -600, Hotspot	A	▲	BRAP p.V600E
Immune Checkpoint Inhibitor	Yes	Signature	A	▲	High tumor mutational burden
MEK inhibitor (Pani, MEK1 inhibitor, MEK2 inhibitor)	Yes	Hotspot	A	▲	BRAP p.V600E
MEK inhibitor (Pani), MEK1 inhibitor, MEK2 inhibitor, RAF inhibitor (Pani)	Yes	Codon -600, Hotspot	A	▲	BRAP p.V600E
RAF inhibitor (Pani)	Yes	Codon -600, Hotspot	A	▲	BRAP p.V600E
Act inhibitor (Pani)	No	Deletion	A	▲	PTEN partial loss
MEK1 inhibitor, MEK2 inhibitor	No	Hotspot	A	▲	BRAP p.V600E
PROCR inhibitor	No	Deletion	A	▲	PTEN partial loss

Potentially relevant clinical studies (NL)

NCT ID	TRIAL	DRUG	MATCH	GENOMIC EVENT
NCT0268687	MR-3475-158KEYNOTE-158	Pembrolizumab	Signature	High tumor mutational burden
NCT02670705	Sub	Atezolizumab, atezolizumab + FOLFIRI	Signature	High tumor mutational burden
NCT0475485	BREFFLY-1	Tovorafenib	Activation	BRAP p.V600E

The Genomemon Clinical Knowledgebase (CKB) is used to generate genomic events with high level evidence. Only evidence of level A (FDA approved therapy and/or guidelines, level B (see clinical trials), and/or level C (early clinical trials) are reported. Evidence level of level D (case reports and preclinical evidence) are not reported. The response symbol ▲ means that the evidence is responsive. The resistant symbol ▼ means that the evidence is resistant. The alteration P (potential) after the response symbol indicates the evidence is predicted responsive/resistant (meaning, the evidence data are limited but a potential responsiveness is suggested). More details about CKB can be found in our Glossary Of Terms.

The Genomemon Clinical Knowledgebase (CKB) database is used to generate genomic events for potential clinical study eligibility. Please note that clinical study eligibility depends on multiple patient and tumor characteristics of which only the genomic events are considered in this report. Therefore, although a clinical study is shown, the patient may not be eligible. When the information is available, clinical studies are filtered on gender.

If the high level evidence or clinical study matching is based on a mutation, but this is not a hotspot (see table Tumor observed variants under Genomic events), that high level evidence or clinical study should be interpreted with extra caution.

If the high level evidence or clinical study matching is based on an amplification, the high level evidence or clinical study that corresponds with 'overexpression' of the gene is also matched. The same rule applies for deletions and underexpression.

For MEK1 genes only: if the high level evidence or clinical study matching is based on an inactivation or deletion, the high level evidence or clinical study that corresponds with 'absence of protein expression' of that gene is also matched.

Using WGS data of the tumor and the reference sample, the molecular tumor cell purity and the average tumor ploidy are estimated.

Tumor specific variants are reported for more than 460 cancer related genes. Only non-synonymous variants are reported and are sorted according to the oncogenic driver likelihood (high, medium and low). Gene coding and protein annotation (VARIANT) is based on the observed chromosomal variants (POSITION) is based on the canonical transcript of the gene and, for certain genes, based on the clinical most relevant transcript. A complete list of the transcripts used can be found in <https://www.oncoact.nl/specheerOncoActWGS>. The READ DEPTH provides the 'raw' sequencing read count of the variant and the total reads observed on the chromosomal position. The tumor variant allele frequency (TVAF) and the gene copy number for all variants have been corrected based on the tumor purity to only represent a tumor specific value. The BIALLELIC column provides information on whether the observed variant is detected in both alleles (bi-allelic) or whether a wildtype allele is still present. A HOTSPOT status highlights the clinical importance of this variant and is provided based on information available from different knowledgebases including CIVIC, DoCM and CGI.

Overview of all observed tumor specific gene disruptions due to structural variants. For each disruption, the disrupted canonical transcript range is shown, as well as the type of disruption (deletions (DEL), inversions (INV), duplications (DUP) and single breaks (BND)) and the number of disrupted and undistrupted allele copies.

Hartwig Medical OncoAct

HOSPITAL PATIENT ID
reporting_normal

REPORT DATE
30-Apr-2023

Genomic events (1/2)

Tumor purity & ploidy

Tumor purity: 99%

Average tumor ploidy: 3.1

Tumor observed variants

GENE	POSITION	VARIANT	READ DEPTH	COPIES	TVAF	BIALLELIC	HOTSPOT	DRIVER
BRAP	7:140453136	c.1799T>A (p.Val600Glu)	150 / 221	6.0	68%	No	Yes	High
CDKN2A	9:219711153	c.246_247delC (p.Gly59fs)	99 / 99	2.0	100%	Yes	Near	High
CDKN2A	9:15	c.203_204delC (p.Asn68fs)	99 / 99	2.0	100%	Yes	Near	High
TERT	5:129528	c.-105_124delC (p.Trp43fs)	56 / 65	2.0	87%	Yes	Yes	High
SF3B1	2:185628779	c.2193C>T (p.Pro714Leu)	74 / 111	3.0	67%	No	No	Low
TP53	17:319604330	c.149T>G (p.Leu59fs)	47 / 112	4.0	42%	No	No	Low

Variant annotation is by default based on the canonical transcript. In case another transcript is more commonly used in routine practice, this annotation is also provided (bottom annotation).

Tumor observed gains & losses

CHROMOSOME	REGION	GENE	TYPE	MIN COPIES	MAX COPIES	CHROMOSOME WITH COPIES
10	q23.31	PTEN	partial loss	0	2	2

Tumor observed gene fusions

None

Tumor observed homozygous disruptions

Complex loss of wild type alleles

None

Tumor observed gene disruptions

LOCATION	GENE	DISRUPTED RANGE	TYPE	CLUSTER ID	DISRUPTED COPIES	UNDISRUPTED COPIES
10q23.31	PTEN	Intron 5 > Intron 6	DEL	68	2.0	0.0

Tumor specific copy number alterations are listed here, including gene copy-gains (amplification) and complete losses.

Gene copy gains are reported if the complete gene (full gain) or only part of the gene (partial gain) shows an increase in copy number. A distinction is made between a sufficiently high enough (defined as higher than 3x the tumor ploidy).

For gene copy losses, only tumor-specific complete losses are reported (0 copies). A distinction is made between a partial loss (only part of the gene has 0 copies) and a full loss (the complete gene has 0 copies) of the gene.

The detected gene fusions that are predicted to result in a viable fusion product are listed here. Information about the fusion partners include:

- The genetic breakpoints of the genes involved (exon level) and their position (5' or 3') in the fusion
- The phasing of the genes ('inframe' or 'out of frame' splicing, which is required for an inframe fusion product)
- The calculated copies of the fusion in the tumor
- The driver likelihood of the gene fusion, with a high-driver status for all known fusions.

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
Using WGS data, the HR status of the tumor can be accurately predicted by the CHORD classifier tool based on specific single nucleotide variants (SNV), insertions and deletions (indels), and structural variant (SV) types. A score higher than 0.5 indicates HR deficiency caused by complete (bi-allelic) inactivation of BRCA1/2 or possibly other genes in the HR pathway (e.g. RAD51C, PALB2). More details are described in Nguyen *et al.* Nature Communications, 2020.

Tumor genomic profiles (1/2)

Homologous recombination status

Proficient 0.00

Microsatellite status

Stable 0.1

Tumors with a microsatellite stability score lower than 4 are considered microsatellite stable (MSS) and tumors with a score larger than 4 are considered microsatellite unstable (MSI). The WGS-based MSI readout has been validated against the routine MSI-POR assay and immunohistochemistry status of proteins involved in the mismatch repair (MMR) pathway.

The tumor mutational burden is reported as:

- The mutational load (ML), which is defined by the total number of somatic missense variants across the whole genome of the tumor.
- Tumor mutational burden (TMB) score, which is calculated by the number of all somatic variants per genome Mb.

Although closely related, differences between both metrics exist. For TMB, tumors with a score >10 are considered to have a high mutational burden, which has clinical significance for possible treatment with immunotherapy.

Tumor mutational burden

High 13.7

Tumor mutational load

165

Manual de usuario de OncoAct

HMF-IVDD-430 V2.5


Hartwig Medical OncoAct

Tumor genomic profiles (2/2)

Molecular tissue of origin prediction

The title shows the conclusion of the prediction of the molecular tissue of origin. If none of the competing predictions has a likelihood $\geq 80\%$, no reliable conclusion can be drawn (results inconclusive).

The left plot shows the likelihoods (percentages) for all the origin types analyzed by the molecular tissue of origin prediction tool. Only when the likelihood is $\geq 80\%$ (a peak in the green outer band of the plot), a reliable prediction (with $\geq 80\%$ accuracy) can be drawn. Lower likelihoods ($< 80\%$) suggest there is a similarity with that tissue of origin, but this is less strong and there is lower confidence.

The right plot shows the breakdown of the strongest predicted likelihood(s) into the contribution of the 11 DNA types related to those used in Cosmic signatures: C1 driver landscape and passenger characteristics (e.g. tumor type observed sharing), and C2 somatic mutation pattern (mutation distribution across the genome).

6/11 HMF-FOR-080 v0.0.7

Hartwig Medical OncoAct

CIRCOS plot

The outer first circle shows the chromosomes. The darker shaded areas represent large gaps in the human reference genome: i.e. regions of centromeres, heterochromatin & missing short arms.

The second circle shows all tumor observed variants (not, exon, intron and intergenic regions) and are divided into an outer ring of single nucleotide polymorphism (SNP) allele frequencies and an inner ring of short insertion/deletion (indel) locations. Variant allele frequencies have been converted to tumor purity and scale from 0 to 100%. Each dot represents a single variant and are colored according to the type of base change (e.g. C to T/G to A in red) and are in concordance with the coding used in Alexandrov et al. (2013) Nature paper that describes the use of mutational signatures. INDELs are colored yellow and red for insertions and deletions respectively.

The third circle shows all observed tumor purity adjusted copy number changes, including both focal and chromosomal events. Copy number losses are indicated in red, green shows regions of copy number gain. The scale ranges from 1 (complete loss) to 6 (high level gains). If the allelic copy number is > 6 it is shown as 6 with a green dot on the diagram.

The fourth circle represents the observed minor allele copy numbers across the chromosome. The region of the chart is from 0 to 3. The expected normal minor allele copy number is 1, and anything below 1 is shown as a bar and represents a CNV event (deletion). Minor allele copy numbers above 1 indicate amplification events of both A and B alleles at the indicated locations (blue).

The innermost circle displays the observed structural variants within or between the chromosomes. Translocations are indicated in blue, deletions in red, inversions in yellow, tandem duplications in green and insertions in black.

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A CIRCOS plot visualizes the position, size and orientation of all tumor-specific genomic alterations. These alterations include single nucleotide variants, insertions/deletions, copy numbers changes, translocations and other structural variants.

Details on how to interpret a CIRCOS plot are described on the bottom of this page.

Manual de usuario de OncoAct

HMF-IVDD-430 V2.5

[illegible]

Each OncoAct report ends with more Sample details and Disclaimers (page 10) and the signature of the Director Hartwig Medical Foundation (page 11).

The Sample details include additional information about the sample processing and report generation. The Disclaimers include the version of the report, the UDI-DI of the OncoAct product and general aspects of the performance of OncoAct (sensitivity).

For feedback and complaints, please contact: qualitysystem@hartwigmedicalfoundation.nl.

For questions regarding the contents of a report, please contact: diagnosticssupport@hartwigmedicalfoundation.nl.

 Hartwig Medical OncoAct	 Hartwig Medical OncoAct	HOSPITAL PATIENT ID reporting_normal REPORT DATE 30-Apr-2025
Sample details & disclaimers (1/2)		
Sample details The samples have been sequenced at Hartwig Medical Foundation, Science Park 408, 10800 Amsterdam. The hospital patient ID is: reporting_normal and the pathology issue ID is: pathologynumber. The results in this report have been obtained between 01-Jan-2023 and 30-Apr-2023. This analysis is performed on the tumor sample as arrived on 01-Jan-2023 with barcode normal. The results stated in this report are based on the tested tumor and reference sample. This equipment is performed according to lab procedures: ISO17 v4.1-SPECT v4.0-PREPRINT v4.0 This report is performed to: studyR, official@hospitalname, hospitalPostCode.hospitalCity <i>Comments: this is a test report based on the COLO829 cell line.</i>	Disclaimers The data of which this report is based is generated from tests that are performed under NEN-EN-ISO/IEC 17025:2017 TESTING L&D accreditation and have passed all internal quality controls. This report is generated using the molecular pipeline version 5.39 and OncoAct reporting engine version 1.6. ID: HRM20230500000APR03-1-1 This report was generated automatically and checked by a trained Clinical Molecular Biologist in Pathology (NMMP). The primary tumor location and primary tumor type have influence on the clinical end-use utility matching. No check is performed to verify the received information. The conclusion of this report is based solely on the results of the whole genome sequencing of the received biospecimen, and the additional primary tumor location and type information received from the hospital. Further interpretation of these results with the patient's clinical context is recommended by a clinician with support of a molecular tumor board. Based on a limited tumor purity of at least 20%, the test has a sensitivity of ~ 95% for detection of tumor observed variants, tumor observed gene and losses, tumor observed gene fusions and tumor observed gene/homologous deletions. Hartwig Medical Foundation is not responsible for the content of all external data sources used to do the analyses and generate this report. Hartwig Medical Foundation is not liable and cannot be held accountable for any inaccuracies, incompleteness or error of any other kind in these data sources, or the external software used to harmonize and curate these data sources. Based on the Dutch Act on Exceptional Medical Treatments (in Dutch: "Wet bijzondere medische verrichtingen") (Storting Hartwig Medical Foundation is not authorized to interpret germline results, but may share data on explicit request of a licensed clinical genetics laboratory that is authorized to provide germline test results. For feedback or complaints please contact: qualitycontrol@hartwigmedicalfoundation.nl. For questions about the contents of this report, please contact: diagnosticssupport@hartwigmedicalfoundation.nl.	— End of report —
	 Edwin Cuppen, Director Hartwig Medical Foundation	